

**Experience Counts**

*With our experience in working laboratories , we understand the lab must be efficient. Whether building a new or modifying a lab, it is more complex. There are many types of work, staffing levels, instrumentation needs and accommodation needs for all who work in the lab space. Staff must participate in layout, configuration of space & necessary instrumentation lists, work flow planning, storage requirements, transports, cleaning, special power needs, lighting, computer & desk & chair needs and area identification. Special requirements for microscopes, biohazards, heat generation & equipment: (fume hoods, glassware washers, storage, walls, tall storage , islands, tables, carts & seating) & Safety is a lot to plan for. We try to ask the right questions, to void errors or omissions While labs are our specialty, each is unique, our experience can provide solutions when a problem interferes with your lab's function or how it will incorporate change.*



**Active SAM Registration  
DUNS# 102851987  
CAGE CODE# 1GU97  
EIN 91-2041685**

**NAICS Codes:  
337127  
Benches Institutional furniture**

**423490  
Laboratory Equipment wholesaler**

**238390  
Lab furniture installations**

**Small Business -(WOSB)**

**Distributor of leading fume hood,  
furniture manufacturers.**

**Contact: Barbara Rodgers  
Phone: 509-796-2599  
Fax: tone activated or  
Email: bdlabcab@ieway.com  
www.laboratoryinteriors.com  
www.lab-specialties.com  
www.ibuildlabs.com**

**Declaration of  
Capabilities:  
Bonded Licensed  
BDCABI\*011CW**

**BD Cabinets, Inc dba  
Laboratory Interiors  
PO Box 335  
Davenport WA 99122  
509-796-2599**



Need stainless steel sloped top tall storage cabinets for a cleanroom or need other lab specials, we can help.

**l a b - s p e c i a l t i e s**

**Safety Cabinets for acids, bases, flammables in stock & custom**



Sliding doors that have self closing feature are required by code in some places. Code and compliance is part of every project or purchase.

Labeling is required. What once was acceptable may not be as newer scientific principles are applied. The best recommended practice for venting flammable or acid cabinets is now a 2" isolated metal for flammables duct to the outside. Venting into other units or hoods is no longer an accepted by all health and safety professionals.

Wood bases can absorb and become saturated & are not preferred for acid or flammable storage unless fully plastic lined and with a pan bottom.

Special base unit for fume hood chemical storage. Acid bases are often now polypropylene.

Lab specials are our specialty.  
[www.lab-specialties.com](http://www.lab-specialties.com)

We offer a plan, specifications or costing review service. To avoid errors or omissions this is available as a flash or comprehensive review service  
Minimum fee for flash plan or specifications review is \$150 a room or \$100 for a specification item.

*Lab Furnishing include:  
Ergonomic mats*

*Configurations of Workstations*

*Laboratory Benches & Casework*

*Balance Tables*

*Storage Cases*

*Fume Hoods*

*Chemical Storage*

*Fixtures & Sinks*

*Emergency showers & Eye wash station*

*Tables, Carts*

*Adjustable Tech Workstations & Desks*

*Glassware washing equipment,*

*Stainless Steel Counters and Casework*

*Epoxy Resin & Trespa Counters & Shelving*

*Lab Stools & Ergonomic Laboratory Chairs*

*Fans for Hoods & IVT (instrument vents)*

**Customers:**

State of WA/Tables **Universities:**

Washington State University– Hoods

University of Missouri Reactor– polypro

University of Washington-Cleanroom

**Wineries:** Hogue Cellar, Columbia Crest, cabinets & epoxy counters

**Chemical Companies Cabinets:**

EkaNobel-

Waste Management

**Water & Waste Water labs D&B**

**Dams:** McNary, John Day, & others

**Fisheries** – Mt. Baker, Leavenworth .

**BD Cabinets, Inc.**

**Barbara Rodgers  
Product Technical Consulting  
509 796-2599**

[www.laboratoryinteriors.com](http://www.laboratoryinteriors.com)  
Email us at [laboratoryinteriors@gmail.com](mailto:laboratoryinteriors@gmail.com)